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Reference	Agomir & Antagomir	Animal	Delivery	Doses	Detection
Jie Ding, <i>et al. Nat. Cell Biol.</i> 2010	hsa-miR-151a-5p antagomir	mice	intrahepatic injection with antagomir-151-5p treated GFP-labeled MHCC-LM3 cells	N/A	5 weeks after injection
Fei Wang, <i>et al. J. Cell. Biochem.</i> 2012	rno-miR-21-5p antagomir	rats	applied with the F-127 pluronic gel	0.24 mg/200 mL 20% F-127 pluronic gel	21 days after injury
Tao Li, <i>et al. J. Urol.</i> 2012	hsa-miR-21-5p antagomir	mice	intratumoral injection	0.2 mg, 3 times per week for 2 weeks.	35 days after the first injection
Minfeng Shu, <i>et al. Mol. Cancer.</i> 2011	hsa-miR-335-5p antagomir	mice	intratumoral injection	0.2 mg, every two days for 2 weeks	3 weeks after the first injection
Yanjie Lu, <i>et al. Circulation.</i> 2010	mmu-miR328-3p antagomir	mice	tail vein injection	80 mg per kg body weight, 1 to 3 consecutive days	0-14 days after the first injection
Hui Li, <i>et al. J. Clin. Invest.</i> 2009	mmu-miR-2861 antagomir	mice	tail vein injection	80 mg per kg body weight, 1 to 3 consecutive days	4 days, 3 weeks and 6 weeks after the first injection
Xiaogang Wang, <i>et al. Nat. Med.</i> 2013	mmu-miR-214-3p antagomir	mice	tail vein injection with bone-targeting delivery system	10 mg per kg body weight, 1 to 3 consecutive days, or every two weeks in 6 weeks	2 weeks or one month after the last injection
Tao Wang, <i>et al. Am. J. Pathol.</i> 2012	mmu-miR-21a-5p antagomir	mice	dermis injection	16 µg per wound for one time	0-18 days after the first injection
Jian-Zhong Hu, <i>et al. J. Neurotrauma.</i> 2013	rno-miR-21-5p antagomir	rats	intrathecal injection	20 nmol/ mL, 1 µL/h for 3 days	4 weeks after injection
Linhui Liang, <i>et al. Hepatology.</i> 2010	hsa-miR-125b-5p antagomir	mice	subcutaneous transplantation with antagomir-125b treated SK-Hep-1 cells	N/A	4 weeks after transplantation
Yanxin Chang, <i>et al. J. Hepatol.</i> 2013	hsa-miR-20a-5p antagomir	mice	subcutaneous or spleen transplantation with antagomir-20a treated GBC-SD cells	200 nM	2 weeks after transplantation
Meijuan Zhou, <i>et al. Carcinogenesis.</i> 2013	hsa-miR-365-3p antagomir	mice	subcutaneous transplantation with antagomir-365 treated A431 cells and intratumoral injection with antagomir-365	100 nM, treated for 24h; 2.5 pmol per tumor mass 3 times per week for two weeks	21 days after treatment
San-Jian Yu, <i>et al. Clin. Cancer Res.</i> 2013	hsa-miR-200a-3p antagomir	mice	tail vein injection with antagomir-200a treated MCF-7 cells	N/A	3 weeks after injection
Mei Wang, <i>et al. Eur. J. Cancer.</i> 2013	hsa-miR-17-5p & hsa-miR-20a-5p antagomir	mice	intratumoral injection	N/A	20 days after injection
Wang-Yu Cai, <i>et al. J. Cell Sci.</i> 2013	mmu-let-7a-5p agomir	mice	intramammary infusion	5 nmol every three days for 4 times	20 days after the first injection
Chuanliang Xu, <i>et al. Mol. Cancer Ther.</i> 2013	hsa-miR-100-5p agomir	mice	intratumoral injection	2 nmol every three days for 4 weeks	5 weeks after injection
Ping Wang, <i>et al. Mol. Cell. Biol.</i> 2013	Mmu-miR-329-3p	mice	intravitreal injection	5 µg	5 days after injection
Yang Sun, <i>et al. Cell Res.</i> 2013	mmu-miR-124-3p agomir	mice	tail vein injection	100 nmol/kg for 3 days	24 hours after the last injection